



Diamantscheiben / Diamond discs

Ausführung / style FLEX

(355)



REF	D190F	D190F
Ø (1/10 mm)	190	190
C (in mm)	0,17	0,25
Körnung / Grit		
Belegung / coating	3,0	3,0
Schaft / shank	HP	HP
	1	1

D220F	D220F
220	220
0,17	0,25
3,0	3,0
HP	HP
1	1




empfohlene Drehzahl Diamantscheiben /
recommended rotation speed diamond discs

HP-Schaft: 20.000 min⁻¹ /
HP shank: 20.000 rpm

Alle Diamantscheiben sind auch unmontiert (lose) lieferbar (5 Stk/Pck).
All diamond discs are also available unmounted (5 pieces/pack).

Ausführung / style SUPERFLEX

(355)






REF	D160S	D160S
Ø (1/10 mm)	160	160
C (in mm)	0,17	0,25
Körnung / Grit		
Belegung / coating	3,0	3,0
Schaft / shank	HP	HP
	1	1

D190S	D190S	D190S
190	190	190
0,14	0,17	0,25
3,0	3,0	3,0
HP	HP	HP
1	1	1

D220S	D220S	D220S
220	220	220
0,14	0,17	0,25
3,0	3,0	3,0
HP	HP	HP
1	1	1

Ausführung / style ULTRAFLEX

(355)



REF	D160U	D160U	D160U
Ø (1/10 mm)	160	160	160
C (in mm)	0,14	0,17	0,25
Körnung / Grit			
Belegung / coating	3,0	3,0	3,0
Schaft / shank	HP	HP	HP
	1	1	1

D190U	D190U	D190U
190	190	190
0,14	0,17	0,25
3,0	3,0	3,0
HP	HP	HP
1	1	1

D220U	D220U	D220U
220	220	220
0,14	0,17	0,25
3,0	3,0	3,0
HP	HP	HP
1	1	1

Diamant-Trennscheiben / Diamond cutting discs (361)

Ausführung / style SUPERFLEX



REF	D080 S 015
Ø (1/10 mm)	80
C (in mm)	0,15
Körnung / Grit	
Belegung / coating	2,5
Schaft / shank	HP
	1

D100 S 015
100
0,15
3,6
HP
1

30.000 min⁻¹

Körnungen / Types of Grit

- Superfein
superfine (15 µm)
- Fein
fine (30 µm)
- Mittel
medium (45 µm)