


HP-Schraub-Mandrell	RA-Schraub-Mandrell	HP-Schraub-Mandrell	HP-Sandpapier-Mandrell, konisch, vernickelt	HP-Schraub-Mandrell 305/080 (verstärkt)	HP-Spezial-Mandrell für Occlupol 2 mm	HP-Spezial-Mandrell für Occlupol 3 mm
<i>HP-Screw mandrel</i>	<i>RA-Screw mandrel</i>	<i>HP-Screw mandrel</i>	<i>HP-sandpaper mandrel, conical, nickel-plated</i>	<i>HP-Screw mandrel 305/080 (strengthened)</i>	<i>HP-Special mandrel for Occlupol 2 mm</i>	<i>HP-Special mandrel for Occlupol 3 mm</i>



REF	PA4001	PA4004	PA4005	PA4007	PA4013	PA4029	PA2007	PA3007
	5/1	5/1	5/1	5/1	5/1	5/1	5/1	5/1

REINIGUNGSSTEIN

zum Reinigen von Gummipolierern und zur Reinigung von Diamantinstrumenten

CLEANING STONE

for cleaning of rubber polishers and diamond instruments




REF	PA9824
	1

DIAMANTIERTER ABRICHTSTEIN

- zum Zentrieren und Zurichten von Gummipolierern und Steinen
- verschiedene Winkel und Radien sowie ebene Flächen

DIAMOND DRESSING STICK

- for centering and shaping of rubber polishers and stones
- different angles and radii as well as flat structures



REF	PA4060
	1