


Peeso


Drehzahl 800 - 1.200 min⁻¹



Schaftkodierung / Ringe		2	3	3
REF 6 Stk./Pck.	WPE 1830...	192	193	194
Ø in 1/100 mm		090	110	130
L in mm		19	19	19

Vorböhrer


Drehzahl 500 - 1.000 min⁻¹



Farbkodierung		1	1	2	2	3	3	4	4
REF	3 Stk / Pck	6320	6321	6322	6323	6324	6325	6326	6327
Bohrerkopf Ø	Ø in mm	1,2	1,2	1,4	1,4	1,4	1,4	1,65	1,65
Bohrerlänge	L in mm	9	16	10	16	10	16	12	16


Planfräser

Drehzahl 2.000 min⁻¹



Farbkodierung		1	2	3	4
REF	3 Stk / Pck	6330	6331	6332	6333
Zapfen Ø	D in mm	1,1	1,3	1,3	1,55
Zapfenlänge	L in mm	4	4	4	4
Kopf Ø	DH in mm	3	3,5	4	4,5

Steckschlüssel für HeadMaster konisch



Farbkodierung	1	2	3	4
REF 1 Stk	6350	6351	6352	6353